



3-JOINT-ARM FRICTION



Product description

3-joint-arm friction with knob for use with semirigid endoscopes:

- fixes the ocular with which is connected to the optic body and the light cable
- can be fixed in every position and serves as a "third" hand
- keeps all accessories in the unsterile area

Order No.

PD-AC-0092

Material

Aluminium

Weight

1,63 kg

Utilisation environment

Operating conditions at room temperature:

15 – 35 °C (59 – 95 °F)

Rel. humidity: 30 – 60 %

Air pressure: 760 – 1060 mbar

Cleansing and disinfection / sterilisation

Wipe disinfection

Certification

CE: This device satisfies the requirements of Council Directive 93/42/EEC as regards medical devices.

This device is a Class I product.