



POLYDIAGNOST

MICROPERC[®] SET



Product description

The MICROPERC[®] set for percutaneous micro endoscopic diagnosis and therapy of kidney stones is exclusively intended for use with semirigid optics (PD-PN-0080) from POLYDIAGNOST GmbH.


The valuable MicroPerc[®] optic comes with the proven PolyScope lens system with its field of view of 120 degrees, reflecting peripheral and direct human vision. MicroPerc[®] can be connected to all conventional endoscopy towers making it easy and economical to expand your toolbox for kidney stone management.

Set consists of:

Introducer needle
8 F shaft
4.8 shaft
Dilator for the 8 F shaft
3-port-adapter for 4.8 F and 8 F shaft
Cap
Tuohy borst adapter
3-way stop cock
Perfusion tube
Lunderquist guidewire



POLYDIAGNOST

Order-No.	PD-PN-1008
Material	Polycarbonate Polypropylene Stainless steel
Weight	0,092 kg
Utilisation environment	Operating conditions at room temperature: 15 – 35 °C (59 – 95 °F) Rel. humidity: 30 – 60 % Air pressure: 760 – 1060 mbar
Cleansing and disinfection / sterilisation	Single use
Packaging	Sterile packed
Certification	CE: This device satisfies the requirements of Council Directive 93/42/EEC as regards medical devices. This device is a Class IIa product.
	PMT GmbH Im Gewerbegebiet 11 D – 66709 Weiskirchen

POLYDIAGNOST GmbH
Am Söldnermoos 17
85399 Hallbergmoos

www.polydiagnost.com